
S1-E / S1-A

**Sendoline<sup>®</sup>**

English / Svenska / Norsk / Dansk


Instructions for use  
Bruksanvisning

Bruksanvisning  
Brugsanvisning

S1-E / S1-A

**Sendoline<sup>®</sup>**

English

Instructions for use

# Content

<b>1</b>	<b>Before you begin</b> .....	<b>5</b>
	1.1 Structure of the document .....	6
<b>2</b>	<b>Safety instructions</b> .....	<b>9</b>
<b>3</b>	<b>Technical description</b> .....	<b>11</b>
	3.1 Task .....	11
	3.2 Set up .....	12
	3.3 Overview of files .....	13
	3.4 Technical Data .....	15
	3.5 Product labeling .....	18
<b>4</b>	<b>Preparation</b> .....	<b>20</b>
	4.1 Initial start-up and longer breaks in use .....	20

<b>5</b>	<b>Operation .....</b>	<b>21</b>
5.1	Replacing the instrument .....	22
5.2	Attaching and removing a file .....	23
5.3	Attaching the Apex clip .....	25
5.4	Connecting the instrument to an endometer.....	26
5.5	Recommendation for the treatment procedure .....	28
5.6	Fitting and removing an insulating sleeve.....	30
<b>6</b>	<b>Postprocessing .....</b>	<b>32</b>
6.1	After each treatment session .....	32
6.2	At the end of the work day .....	33
<b>7</b>	<b>Reprocessing .....</b>	<b>34</b>
7.1	Conduct pre-disinfection .....	34
7.2	Automated cleaning and disinfection .....	35

	7.3	Manual cleaning and disinfection.....	38
	7.4	Manual maintenance .....	40
	7.5	Sterilizing .....	44
8		Spare parts and consumables .....	46
9		Disposal .....	47

## 1 Before you begin ...

S1-E / S1-A complies with the regulations taking into account the state of the art. S1-E / S1-A complies with the ISO 14457 standard.

1. Read the operating instructions prior to using S1-E / S1-A.
2. Use S1-E / S1-A only for the applications described in the operating instructions.
3. Observe the hygiene standards, occupational safety regulations, and accident prevention measures applicable to the S1-E / S1-A in your country.

S1-E / S1-A is used to hold and drive burs for the purposes of rotary processing and is intended for the following uses in dentistry:

- Endodontics
- Endometry

none

### Intended use

### Contraindications

## Target group

This product is intended only for use by trained dental personnel in dental practices and laboratories.

### 1.1 Structure of the document

#### 1.1.1 Labeling of information

## Warnings

➤ To prevent injuries, please observe the warnings.

Warnings are labeled as follows:

⚠ **DANGER!** indicates a danger **leading** to death or serious injury if not avoided.

⚠ **WARNING!** indicates a danger that **may lead** to death or serious injury if not avoided.

⚠ **CAUTION!** indicates a danger that **may lead** to injury if not avoided.

## Instructions for use

➤ To prevent material damage and additional expenses, please observe all instructions for use.

Instructions for use are labeled as follows:

**NOTE!** indicates measures for the prevention of material damage.

**IMPORTANT:** indicates information on the avoidance of additional expenses and other important information.

**Tip:** indicates information for facilitating work.

### 1.1.2 Formats and symbols

The formats and symbols used in this document have the following meaning:

<ul style="list-style-type: none"> <li>✓ Prerequisite</li> <li>1. First action step</li> <li>2. Second action step</li> <li>or</li> <li>➤ Alternative action</li> <li>↩ Result</li> <li>➤ Individual action step</li> </ul>	<p>Requests you to do something.</p>
<p>Use of formats and symbols [→ 8].</p>	<p>Identifies a reference to another text passage and specifies its page number.</p>
<ul style="list-style-type: none"> <li>• List</li> </ul>	<p>Identifies a list.</p>

## Obligations of the user

## Preventing the spread of infections and cross contamination

## 2 Safety instructions

- Use only fault-free materials that do **not** deviate from the specified data [→ 15].
- Protect yourself, patients, and others from danger. To do this, observe the safety information.
- Comply with the Intended use of the equipment.
- You should always keep these operating instructions within reach for further reference.

Prevent the spread of infections and cross contamination between patients, users, and third parties. Sterilize equipment after each patient.

Take the appropriate hygiene measures, e.g. wear protective gloves.

## **Instrument head/body overheating**

If the instrument is defective, the area around the instrument head or the instrument body may heat up, thus creating a risk of burning the hands or the patient's oral mucosa.

## **Malfunction or damage**

Discontinue use immediately in case of malfunction, unusual or different sounds or damage. Damaged instruments may cause injury. Notify the dental depot or the manufacturer.

## **Repair**

Do **not** repair the instrument yourself.

## **Spare and accessory parts**

Use only original parts produced by the manufacturer.

## **Operating conditions:**


Temperature: +10 °C – +30 °C.

If you have any questions, please contact your dental depot or the manufacturer.

## 3 Technical description


### 3.1 Task

The instrument serves for transmitting the driving power of the electric motor or air motor (handpiece coupling according to ISO 3964) to the preparation tool and converts a rotational movement into a rotating movement with alternating direction.

## 3.2 Set up

A	Instrument head
B	Push button
C	Opening of chuck system
D	Joint for clip

## Accessories

A	Insulating sleeve
B	Clip

### 3.3 Overview of files

Following files are to use for the application with S1-E / S1-A:

Designation	ISO/Taper
Sendoline - S1 Plus Standard	25/06
Sendoline - S1 Plus Small	15/06

<b>Designation</b>	<b>ISO/Taper</b>
Sendoline - S1 Plus Large	40/04
Sendoline – S1 File Standard 25/06	25/06
Sendoline – S1 File 15/06	15/06
Sendoline - S5 Rotary File #1	30/08
Sendoline - S5 Rotary File #2	30/06
Sendoline - S5 Rotary File #3	30/04
Sendoline - S5 Rotary File #4	25/04
Sendoline - S5 Rotary File #5	20/04
Sendoline - Big Apical File	35/04
Sendoline - Big Apical File	40/04

Designation	ISO/Taper
Sendoline - Big Apical File	45/04
Sendoline - Big Apical File	50/04

### 3.4 Technical Data

	S1-E	S1-A
Gear ratio	~ 150:1	~ 82:1
Maximum motor speed in rpm	~ 40 000	~ 22 000
Mean speed at maximum motor speed in rpm	~ 270	~ 270
Motion	forward and backward rotating	forward and backward rotating
Clockwise rotation angle	~ 180°	~ 180°

	<b>S1-E</b>	<b>S1-A</b>
Counterclockwise rotation angle	~ 30°	~ 30°
Periodicity of the motion in Hz	~ 10.7	~ 10.7
Clamping system	CA	CA
Internal cooling media	-	-
Light function	-	-
Apex measurement function	x	x
Clips available for standard endometer connections	∅ 2mm	∅ 2mm
Instrument coupling	ISO 3964	ISO 3964

**File**

	<b>S1-E / S1-A</b>
Shank diameter in mm	2.334 - 2.350
Maximum total length in mm	42
Maximum working diameter in mm (ISO 2157)	2.1
Gripping length in mm	$\geq 11$
Standard	ISO 1797-1

### 3.5 Product labeling

Can be thermally  
disinfected

Sterilizable in a  
steam sterilizer  
(autoclave) at  
the temperature  
specified

Data matrix code (here: an example)

Content of the data matrix code:

+ESEN6542471/\$\$+5192/16D20170309E

A

B

C

D

- A Manufacturer code (here: ESEN)
- B REF number (here: 6542471)
- C Serial no. (here: 5192)
- D Manufacture date (YYYYMMDD)

## **4 Preparation**

### **4.1 Initial start-up and longer breaks in use**

- Sterilize the instrument and accessories prior to startup.
- Clean and maintain the instrument after longer breaks in use.

## 5 Operation

**NOTE!** Use only faultless files in order to prevent fatigue fractures of the files as far as possible.

**⚠ CAUTION!** Use the S1-A only with an air-driven motor at a maximum speed of 22,000 rpm.

**⚠ CAUTION!** Never pull the patient's cheek back with the contra-angle handpiece! This would actuate the pushbutton, thus creating a risk of burning the patient's oral mucosa.

**⚠ CAUTION!** Burning hazard at push button! Press the push button only when the instrument has come to a stop.

**IMPORTANT:** For application details and operating data of the different preparation tools, please refer to the information provided by the manufacturer.

## 5.1 Replacing the instrument


**⚠ CAUTION!** The instrument should only be fitted or removed when the motor is at standstill.

### Attaching the S1-E / S1-A

- ✓ The motor is at a standstill.
- Insert the instrument until it snaps into place.

### Removing the S1-E / S1-A

- ✓ The motor is at a standstill.
- Detach the instrument. Do not pull on the supply hose while doing this.

## 5.2 Attaching and removing a file


**NOTE!** Use only faultless files in order to prevent fatigue fractures of the files as far as possible.

**IMPORTANT:** Check the push button to make sure it moves freely!

**IMPORTANT:** For safety reasons, check the clamping system of contra-angle handpiece bur instruments before every use.

**NOTE!** Only use files, which are approved by Sendoline.

**IMPORTANT:** The file is designed for reciprocal use.

### Attaching a file


- ✓ The motor has come to a stop.
- 1. Insert the file by gently rotating it until it snaps into place. In doing so, do **not** press the pushbutton.
- 2. Then pull and turn the file to check that it is firmly seated.

### Removing a file

- ✓ The file has come to a stop.
- Press the push button and detach the file.

### 5.3 Attaching the Apex clip

1. Attach the APEX clip to S1-E / S1-A.
2. Connect an endometer [→ 26].


## 5.4 Connecting the instrument to an endometer

The S1-E / S1-A makes it possible to connect a measuring cable for an endometer by means of a clip in order to check the exact position of the file in the root canal.

**NOTE!** Use only endometers, which are approved in the respective country.

1. Please observe the operating instructions pertaining to the endometer being used.
2. Attach an Apex clip [→ 25].
3. Connect the measuring cable of the endometer to the clip's contact pin.
4. Perform the measurement in accordance with the equipment manufacturer's specifications.

Based on in vitro studies, we recommend defining the maximum preparation depth at the coronal end of the apex interval on the display of the endometer used.


## Calibration

- ✓ The manufacturer of the endometer stipulates that a calibration must be performed.
- 1. Insert a file into the contra-angle handpiece [→ 23].
- 2. Connect the cable of the file clamp (without the actual file clamp) onto the clip of the instrument.
- 3. Perform the calibration in accordance with the equipment manufacturer's specifications.

## 5.5 Recommendation for the treatment procedure


**NOTE!** Use only faultless files in order to prevent fatigue fractures of the files as far as possible.

- The different file systems require different working techniques. For this reason, please always observe the information provided by the file manufacturer.
- Work without exerting any pressure on the file.
- Prior to the first treatment, we recommend performing conditioning tests on resin blocks or extracted teeth.
- Due to the special characteristics of the nickel-titanium material, the crown down technique should be preferred when preparing the root canal.
- The ApexLocator function is used to support the root canal preparation process in combination with a treatment center with the apex measurement function. You must always take at least one X-ray to determine the preparation depth. A clinical

evaluation including knowledge of the anatomy of root canals is important for interpreting the results.

- When performing the apex measurement, wear appropriate insulated gloves. We recommend using a cofferdam for treatment and putting an insulating sleeve on the contra-angle handpiece. This will prevent inaccurate measurements caused by undesired leakage currents. During the measurement, the instrument must not come into contact with the patient's mucosa, metallic tooth restorations, or the mucosal electrode.
- If the file stops moving, remove the file by gently pulling it toward the coronal end. If this is not possible, switch the motor to the opposite rotation direction and pull the file out of the canal.  
**IMPORTANT:** Check the file afterwards for damage or deformation (untwisting) and replace it if necessary.

## 5.6 Fitting and removing an insulating sleeve

**⚠ WARNING!** The insulating sleeve is not sterile! To avoid cross-contamination sterilize the insulating sleeve before use [→ 44].

**NOTE!** The insulating sleeve is intended for single use. Dispose of the insulating sleeve after each patient.

### Fitting an insulating sleeve

- ✓ You have sterilized the insulating sleeve.
- ✓ The motor is at a standstill.
- ✓ A file is **not** clamped in place.

1. Pull the insulating sleeve fully over the instrument by gently turning it from one side to the other.
2. Insert the file [→ 23].
3. Attach the contra-angle handpiece to the motor.

#### **Removing an insulating sleeve**

- ✓ The motor is at a standstill.
1. Pull the contra-angle handpiece off of the motor.
  2. Remove the file.
  3. Remove the insulating sleeve.
  4. Dispose of the insulating sleeve [→ 47].

## 6 Postprocessing

### 6.1 After each treatment session

**NOTE!** Condition immediately, or at the latest, one hour after treatment.

**NOTE!** Only use W&H Service Oil F1, MD-400.

- ✓ Wear appropriate protective clothing.
- 1. Remove the file with tweezers.
- 2. Predisinfect directly at the treatment center [→ 34].
- 3. Detach the instrument from the motor.
- 4. Transport the instrument to the hygiene room in a suitable transport container.
- 5. Perform automatic reprocessing [→ 35]. Manual reprocessing [→ 38] is possible in exceptional cases if the national/local regulations are followed.

6. Apply spray to the instrument [→ 40].
7. Sterilize the instrument [→ 44].

## 6.2 At the end of the work day

- Apply spray to the instrument [→ 40].

**NOTE!** Do not leave any instruments on the motor overnight, in order to prevent oil from leaking into the electric motor. Never lubricate the electric motor.

## 7 Reprocessing


### 7.1 Conduct pre-disinfection

- ✓ Wear appropriate protective clothing.
- ✓ All disinfectants must be approved in your country and have proven bactericidal, fungicidal and virucidal properties. Use only disinfectants with **no** protein-fixing properties.
  1. Wipe the surface with a disinfectant cloth.
  2. Wipe the disinfectant off with a cloth.
- ✎ For further reprocessing, the instrument should be dry and free of residue.

Please observe the manufacturer's instructions for using the instrument disinfectant.

## 7.2 Automated cleaning and disinfection ...

### 7.2.1 ... with a DAC Universal

We recommend using **Dentsply Sirona DAC Universal** for automated cleaning, disinfection, and care.

For further details, refer to the instructions for use supplied with the unit.

- ✓ The pushbutton chuck is lubricated manually.
- ✓ The instrument is conditioned using the DAC Universal.
- 1. Check whether the instrument is clean after reprocessing under good lighting (min. 500 lux) and color rendering index (min. 80 Ra).
- 2. If they are still dirty, repeat the process.
  - ✎ For further reprocessing, the instrument should be dry and free of residue.

3. Pack the instrument in packaging material suitable for sterilization and storage. e.g. paper/plastic composite packaging according to ISO 11607.
4. Perform sterilization [→ 44].

### 7.2.2 ... with cleaning and disinfection equipment

The instrument can also be cleaned and disinfected in a suitable piece of cleaning and disinfection equipment.

The cleaning and disinfection equipment used must be approved by its manufacturer for the cleaning and disinfection of dental instruments and comply with ISO 15883-1/-2 (e.g., 95°C (203°F) and 10 min. holding time).

For further details, refer to the instructions for use supplied with the unit.

- ✓ The instrument is conditioned with a cleaning and disinfection device.

1. Check whether the instrument is clean after reprocessing under good lighting (min. 500 lux) and color rendering index (min. 80 Ra).
2. If they are still dirty, repeat the process.
  - ✎ For further reprocessing, the instrument should be dry and free of residue.
3. Blow the instrument out with max. 3 bar.
4. Lubricate mechanical parts manually [→ 40].
5. Maintain the push button chuck manually.
6. Pack the instrument in packaging material suitable for sterilization and storage. e.g. paper/plastic composite packaging according to ISO 11607.
7. Perform sterilization [→ 44].

### 7.3 Manual cleaning and disinfection

**IMPORTANT:** Manual reprocessing is possible in exceptional cases if the national/local regulations are followed. The national/local regulations are to be checked before.

**NOTE!** Condition immediately, or at the latest, one hour after treatment.


**NOTE! Never** clean in an ultrasonic bath!

**NOTE! Never** immerse in disinfectants!

**IMPORTANT:** Use a soft, clean, and disinfected brush for cleaning.

- ✓ Wear appropriate protective clothing.
- ✓ All disinfectants must be approved in your country and have proven bactericidal, fungicidal and virucidal properties. Use only disinfectants with **no** protein-fixing properties.

1. Brush the instrument under running water ( $< 38\text{ }^{\circ}\text{C}$ ,  $< 100\text{ }^{\circ}\text{F}$ , at least drinking water quality) and good lighting (minimum 500 lux) and color rendering index (minimum 80 Ra) until no more dirt can be seen, for at least 10 seconds.
2. Flush the drive channels with W&H Service Oil F1, MD-400.
3. Conduct thermal disinfection or unwrapped steam sterilization [ $\rightarrow$  35].
4. Lubricate mechanical parts manually [ $\rightarrow$  40].
5. Maintain the push button chuck manually.
6. Pack the instrument in packaging material suitable for sterilization and storage. e.g. paper/plastic composite packaging.
7. Perform sterilization [ $\rightarrow$  44].

## 7.4 Manual maintenance

### 7.4.1 Maintenance of mechanical parts


#### Intervals


- Prior to each sterilization
- After every thermal disinfection without integrated maintenance

#### Required accessories

- |   |                            |
|---|----------------------------|
| A | Spray adapter              |
| B | W&H Service Oil F1, MD-400 |

**NOTE!** Only use W&H Service Oil F1, MD-400.

## Procedure

### With W&H Service Oil F1, MD-400

- ✓ The spray adapter is disinfected.
- 1. Fit the spray cap with the nozzle onto the spray can.
- 2. Shake the spray can before each application.
- 3. Insert the nozzle into the Instrument.
- 4. Hold the instrument firmly.
- 5. Spray the instrument for about 1 second.  
**IMPORTANT:** Hold the spray can upright.
- 6. Wipe any spray that comes out with a disinfection cloth.
- 7. Repeat the process until the spray that runs out is clear.

### With W&H Assistina

- Observe the instructions for use of the manufacturer.


## 7.4.2 Care of the push button chuck


### Interval

- At least once a week

### Only with W&H Service Oil F1, MD-400


- ✓ The instrument is clean and disinfected.
- 1. Fit the spray cap on the spray can.
- 2. Hold the instrument firmly.
- 3. Press the tip of the spray nozzle firmly into the chuck system.
- 4. Spray into the instrument for about 1 second.  
**IMPORTANT:** Hold the spray can upright.
- 5. Wipe any spray that comes out with a disinfection cloth.

### Only with W&H Assistina

- ✓ The instrument is clean and disinfected.
- 1. Fit the adapter onto the maintenance attachment.
- 2. Direct the spray exit nozzle of the adapter downwards.
- 3. Press the bur side of the instrument onto the adapter from below.
- 4. Start the Assistina.
- 5. Press the instrument onto the adaptor for about 10 seconds.
- 6. Remove the instrument from the adapter.
- 7. Wipe any spray that comes out with a disinfection cloth.
- 8. Close the cover of the Assistina and allow it to run for the remainder of the maintenance cycle (approx. 25 seconds).

## 7.5 Sterilizing

- ✓ The instrument is clean and disinfected.
- ✓ The instrument is lubricated.
- ✓ The instrument can be sterilized in packaging suitable for sterilization and storage: paper/plastic composite packaging according to ISO 11607.
- Sterilize the instrument in the steam sterilizer with saturated water vapor.

Overpressure: 2.04 bar (29.59 psi)

Temperature: 134 °C (274 °F)

Holding time: 3 min.

Steam sterilizers that meet the requirements of EN 13060, class B or S, and are also suitable for the sterilization of straight/contra-angle handpieces are approved.

**NOTE!** Do not exceed 140 °C (284 °F), even during the drying phase.

#### After sterilizing

1. Remove the instrument from the steam sterilizer immediately.  
**⚠ CAUTION!** The instrument is hot. Risk of burns!  
**NOTE!** Do **not** attempt to accelerate the cooling process by immersing the instrument in cold water. This can damage your instrument.
2. Store all instruments so that they are protected from contamination.
3. Sterilize again once the storage period has elapsed.

## 8 Spare parts and consumables

	REF		REF
Service Oil F1, MD-400 (package of 6 oil spray cans, per 400 ml)	10940021		S1 apex clip Ø 2.0 mm 65 83 566
Spray cap with nozzle	02038200		S1 insulating sleeve (for S1-E/S1- A) 65 83 541
W&H spray cap with spray nozzle for chuck system	02036100		W&H Assistina adapter for chuck system 02693000


## 9 Disposal

- According to current information, the product does not contain any substances that are hazardous to the environment.
- Disinfect the product prior to disposal.
- Observe the applicable disposal regulations for your area.

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**Sendoline AB**